

UV-2000S ULTRAVIOLET TRANSMITTANCE ANA

Instantaneous UVA/UVB protection factor values of sunscreen samples



ADVANCED

Labsphere's UV-2000S incorporates the latest component and software technology into an industry proven system architecture, to achieve accurate in-vitro SPF/UVA-Protection Factor analysis of sun care products developed to receive the "very high" sun protection label. Driven by rapidly evolving industry requirements to simplify product labeling and new in-vitro methods to validate product UVA Protection, the UV-2000S is designed to comply with recently approved in-vitro methods, such as COLIPA UVA-PF, Boots Star Rating and the US FDA, as well as several pending global standards/methods. The UV-2000S has replaced Labsphere's UV-1000S as the Industry's choice for not only laboratory in-vitro SPF/UVA analysis, but also production floor quality control.

FAST

The UV-2000S rapidly measures the diffuse transmittance of sunscreen samples in the ultraviolet wavelength region from 250 - 450 nm. Labsphere's Spectralon® integrating sphere incorporates a re-optimized xenon flash lamp to provide exceptional diffuse illumination of the product sample and minimize data integration time. New high performance diode array spectrometers coupled by new, advanced fiber optics are optimized at the system level for low stray light with superior wavelength stability and flash-to-flash repeatability.

FEATURES:

- One touch sample analysis, with results in less than five seconds
- Manual stage for accurate sample positioning and pre and post irradiation
- New Wavelength standard that captures six relevant spectral bands
- Easy-to-use menu driven application software
- Simple instrument performance validation routine ensures accurate, repeatable measurements
- Automatic calculation of SPF, UVA to UVB ratio, critical wavelength, COLIPA Method, Revised Boots Star Rating (pending), and the new FDA testing method

IMPROVED

Many improvements are incorporated in the UV-2000S to realize a new industry de facto standard for in-vitro sun care product analysis. System improvements include new spectrometers, Xenon Flash Lamp, optical coupling fibers, optical head positioning mechanism, sample positioning stage and a new, robust software development platform.

The diode array spectrometers feature stable, custom concave diffractive optics for measurement integrity and repeatability, original holographic diffraction gratings, not replicated gratings, peaked for higher efficiency across the wavelength range, and longer pixel arrays for better pixel wavelength resolution.

Illumination is filtered at the integrating sphere to limit total exposure at the sample and to improve stray light performance. A higher flash rate reduces exposure time, minimizing dark current and maximizing dynamic range.

Use of solarization resistant fibers maintains high system throughput over time. Longer fibers filter high order modes to provide cleaner grating illumination improving stray light performance.

EASY-TO-OPERATE

A built-in report function generates essential information at the click of a button. Reports include necessary information such as date, time, operator name, sample identification, and test parameters. Reports are conveniently viewed on your PC, printed, or exported as text to data processing software for further review and analysis.

POWERFUL APPLICATION SOFTWARE

Developed with .NET Framework®, the license controlled UV-2000S Software features different in-vitro measurement methods for UVA/UVB protection factors of sunscreen including the COLIPA, Boots Star and FDA methods. By licensing the included software, users are able to receive free upgrades on their software for one year as new methods are released to the industry. This easy to use Vista® compatible software facilitates capture/archival/retrieval and export of all data including bare substrate data that may impact UVA-PF due to surface roughness.

UV-2000S application software includes an integrated Performance Validation Routine that allows for on-site validation and re-validation to ensure optimum instrument performance. A set of calibrated standards, including a wavelength standard that captures six relevant spectral bands, is included with each UV Transmittance Analyzer.

Included Model Name	Order Number
UV-2000S	AA-00909-000
UV-2000S Control Software	AS-02755-001
25 HelioScreen HD2 Plates	PP-02097-000
25 HelioScreen HD6 Plates	PP-02101-000

System Properties and Performance

Wavelength Range	250 to 450 nm*
Wavelength Accuracy	±1 nm
Bandwidth (FWHM)	<4 nm
Wavelength Step (Data Interval)	1 nm
Optical Geometry	Hemispherical Illumination/0° viewing (d/O)
Integrating Sphere Geometry	Spectralon®
Integrating Sphere Port Area	< 5%
Sample Exposure Area	0.79 cm ²
Lamp	Xenon Flash Lamp
UV Dose Per Measurement Cycle	< 0.2 J/cm ²
Sample Positioning Stage	Manual Stage
Measurement Range	
Transmittance	0-100%
Absorbance	0 - 2.7 A (Dual Doped PMMA Method)
SPF	1 - 50+
Scan Time	< 5 s
Measurement Methods Supported	
Bare Substrate Analysis and Data Archival	Yes
SPF	Yes
UVA/UVB	Yes
Critical Wavelength	Yes
UVA Protection Factor - COLIPA Method (2007a)**	Yes
UVA Protection - Revised Boots Star Rating (2008)**	Yes
UVA Protection Factor - FDA UV1/UVA**	Yes
UVA Protection Factor - ISO	Roadmap Upgrade
Computer Interface	USB
Min computer requirements	1.6 GHz processor, Windows®XP or Vista SVGA 800 x 600 256MB RAM, 400MB free disk space
Power Requirements	110 - 120/220 - 240 VAC, 60/50 Hz
Operating environment	0° - 50°C, 0% - 70% RH (non-condensing)
Dimensions	
With Stage	11H x 22.6D x 12.3W In (27.9H x 56.6D x 31.2W cm)
Without Stage	11H x 12.6D x 12.3W in (27.9H 32.0D x 31.2W cm)
Optional Accessories	
UV-2000S Starter Kit	AS-02796-000
HelioScreen HD2 Plates	PP-02097-000
HelioScreen HD6 Plates	PP-02101-000
Updated Software Subscription	UV-2000S Software

* All system specifications are based on a wavelength range of 290 to 400 nm.

** Requires a separate solar simulator



LIBERO
Testing Equipment
測試設備

理寶科技有限公司 Libero Technology Limited

香港 HongKong T: (852) 2555 8222 F: (852) 2518 0115
 上海 Shanghai T: 86 (21) 5655 8285 F: 86 (21) 5655 7752
 北京 Beijing T: 86 (10) 6219 4585 F: 86 (10) 6219 4577
 青島 Qingdao T: 86 (532) 8584 8516 F: 86 (532) 8584 8516
 廣州 Guangzhou T: 86 (20) 3928 3292 F: 86 (20) 3928 3290
 www.liberohk.com Email: sales@liberohk.com

